




S35G1 rabbit pAb

Catalog No	YP-Ab-09037
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SLC35G1 C10orf60 TMEM20
Immunogen	Synthesized peptide derived from human S35G1 AA range: 14-64
Specificity	This antibody detects endogenous levels of S35G1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S35G1;SLC35G1 C10orf60 TMEM20
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Translocates from the endoplasmic reticulum to the cell membrane in response to a depletion of intracellular calcium and is detected at punctae corresponding to junctions between the endoplasmic reticulum and the cell membrane.
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Contains 2 DUF6 domains.,
Background	This gene encodes a transmembrane protein which is a member of the drug/ metabolite transporter protein superfamily. The encoded protein may play a role in the regulation of calcium levels inside the cell. [provided by RefSeq, Sep 2016]
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night