




UN5CL rabbit pAb

Catalog No	YP-Ab-09038
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UNC5CL ZUD
Immunogen	Synthesized peptide derived from human UN5CL AA range: 80-130
Specificity	This antibody detects endogenous levels of UN5CL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UN5CL;UNC5CL ZUD
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type III membrane protein . Cytoplasm .
Tissue Specificity	Expressed in pancreas, liver and kidney.
Function	function:Inhibits NF-kappa-B-dependent transcription by impairing NF-kappa-B binding to its targets.,similarity:Belongs to the unc-5 family.,similarity:Contains 1 death domain.,similarity:Contains 1 ZU5 domain.,subunit:Interacts with p65/RELA and NFKB1.,tissue specificity:Expressed in pancreas, liver and kidney.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night