




TM14A rabbit pAb

Catalog No	YP-Ab-09052
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	TMEM14A C6orf73 PTD011
Immunogen	Synthesized peptide derived from human TM14A AA range: 6-56
Specificity	This antibody detects endogenous levels of TM14A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM14A;TMEM14A C6orf73 PTD011
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane .
Tissue Specificity	Expressed at significantly higher levels in ovarian cancer tissues than in normal tissues (at protein level).
Function	similarity:Belongs to the UPF0136 (TMEM14) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).