




MPPE1 rabbit pAb

Catalog No	YP-Ab-09056
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	MPPE1 PGAP5 PP579
Immunogen	Synthesized peptide derived from human MPPE1 AA range: 89-139
Specificity	This antibody detects endogenous levels of MPPE1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PP579; Post-GPI attachment to proteins factor 5; PGAP5; Metallophosphoesterase 1; EC:3.1.-.-; MPPE1; MPPE1 PGAP5 PP579
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum-Golgi intermediate compartment membrane ; Multi-pass membrane protein . Golgi apparatus, cis-Golgi network membrane ; Multi-pass membrane protein . Also localizes to endoplasmic reticulum exit site.
Tissue Specificity	Expressed in brain.
Function	cofactor: Binds 2 divalent metal cations., similarity: Belongs to the metallophosphoesterase superfamily. MPPE1 family., tissue specificity: Expressed in brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night