




TPRKB rabbit pAb

Catalog No	YP-Ab-09058
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TPRKB CGI-121 My019
Immunogen	Synthesized peptide derived from human TPRKB AA range: 122-172
Specificity	This antibody detects endogenous levels of TPRKB at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGI 121; My019; PRPK binding protein; TP53RK binding protein; TPRKB; TPRKB CGI-121 My019
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol . Nucleus .
Tissue Specificity	Widely expressed.
Function	similarity: Belongs to the CGI121/TPRKB family., subunit: Interacts with TP53RK/PRPK., tissue specificity: Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night