




FOXQ1 rabbit pAb

Catalog No	YP-Ab-09061
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FOXQ1 HFH1
Immunogen	Synthesized peptide derived from human FOXQ1 AA range: 281-331
Specificity	This antibody detects endogenous levels of FOXQ1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HNF 3/forkhead like protein 1; HFH-1; HFH1; HFH 1; forkhead box Q1; FOXQ1; FOXQ1 HFH1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed predominantly in the stomach, trachea, bladder and salivary gland.
Function	similarity: Contains 1 fork-head DNA-binding domain., tissue specificity: Expressed predominantly in the stomach, trachea, bladder and salivary gland.,
Background	FOXQ1 is a member of the FOX gene family, which is characterized by a conserved 110-amino acid DNA-binding motif called the forkhead or winged helix domain. FOX genes are involved in embryonic development, cell cycle regulation, tissue-specific gene expression, cell signaling, and tumorigenesis (Bieller et al., 2001 [PubMed 11747606]). [supplied by OMIM, May 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night