




TSN1 rabbit pAb

Catalog No	YP-Ab-09069
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TSPAN1
Immunogen	Synthesized peptide derived from human TSN1 AA range: 179-229
Specificity	This antibody detects endogenous levels of TSN1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BCLF 1; BCLF1; Component 3 of promoter of RISC; RCHF1; recombination hotspot associated factor; recombination hotspot binding protein; reh1; REHF1; TBRBP; Translin; TRSLN; TSN; TSN_HUMAN; TSN1; TSPAN1
Observed Band	
Calculated Molecular Weight	27kD
Cell Pathway	Lysosome membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity: Belongs to the tetraspanin (TM4SF) family.,
Background	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night