




CYYR1 rabbit pAb

Catalog No	YP-Ab-09099
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	CYYR1 C21orf95
Immunogen	Synthesized peptide derived from human CYYR1 AA range: 52-102
Specificity	This antibody detects endogenous levels of CYYR1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	cysteine/tyrosine rich 1;C21orf95;CYYR1;CYYR1 C21orf95
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Widely expressed.
Function	similarity:Belongs to the CYYR1 family.,tissue specificity:Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4° over night