




T178A rabbit pAb

Catalog No	YP-Ab-09101
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM178A TMEM178 PSEC0131 UNQ5926/PRO19820
Immunogen	Synthesized peptide derived from human T178A AA range: 207-257
Specificity	This antibody detects endogenous levels of T178A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T178A; TMEM178A TMEM178 PSEC0131 UNQ5926/PRO19820
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity: Belongs to the TMEM178 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4°over night