




IGJ rabbit pAb

Catalog No	YP-Ab-09104
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	IGJ IGCJ
Immunogen	Synthesized peptide derived from human IGJ AA range: 81-131
Specificity	This antibody detects endogenous levels of IGJ at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Joining chain of multimeric IgA and IgM;JCHAIN;JCH;Immunoglobulin J chain; IGCJ;IGJ;IGJ IGCJ
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function:Serves to link two monomer units of either IgM or IgA. In the case of IgM, the J chain-joined dimer is a nucleating unit for the IgM pentamer, and in the case of IgA it induces larger polymers. It also help to bind these immunoglobulins to secretory component.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4°over night