




IRGQ rabbit pAb

Catalog No	YP-Ab-09120
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	IRGQ IRGQ1 FKSG27
Immunogen	Synthesized peptide derived from human IRGQ AA range: 388-438
Specificity	This antibody detects endogenous levels of IRGQ at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FKSG27;IRGQ;IRGQ1;IRGQ IRGQ1 FKSG27
Observed Band	69kD
Calculated Molecular Weight	69kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night