




UTP15 rabbit pAb

Catalog No	YP-Ab-09126
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	UTP15
Immunogen	Synthesized peptide derived from human UTP15 AA range: 322-372
Specificity	This antibody detects endogenous levels of UTP15 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NET21; UTP15
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus . Found predominantly at the fibrillar center. .
Tissue Specificity	
Function	function: Involved in nucleolar processing of pre-18S ribosomal RNA., similarity: Contains 7 WD repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night