




AN18A rabbit pAb

Catalog No	YP-Ab-09128
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ANKRD18A KIAA2015
Immunogen	Synthesized peptide derived from human AN18A AA range: 323-373
Specificity	This antibody detects endogenous levels of AN18A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AN18A;ANKRD18A KIAA2015
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Contains 5 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night