




TTC18 rabbit pAb

Catalog No	YP-Ab-09141
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TTC18
Immunogen	Synthesized peptide derived from human TTC18 AA range: 963-1013
Specificity	This antibody detects endogenous levels of TTC18 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTC18;TPR repeat protein 18
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium, flagellum . Cytoplasm, cytoskeleton, flagellum basal body . Cell projection, cilium . Cytoplasm, cytoskeleton, cilium axoneme . Present all along the flagellum, with a marked signal at the base of the flagellum. .
Tissue Specificity	Expressed in testis.
Function	similarity:Contains 8 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night