




DB125 rabbit pAb

Catalog No	YP-Ab-09146
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DEFB125 DEFB25
Immunogen	Synthesized peptide derived from human DB125 AA range: 4-54
Specificity	This antibody detects endogenous levels of DB125 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DB125;DEFB125 DEFB25
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function:Has antibacterial activity .,similarity:Belongs to the beta-defensin family.,
Background	Defensins are cysteine-rich cationic polypeptides that are important in the host immunologic response to invading microorganisms. The antimicrobial protein encoded by this gene is secreted and is a member of the beta defensin protein family. Beta defensin genes are found in several clusters throughout the genome, with this gene mapping to a cluster at 20p13. [provided by RefSeq, Nov 2014],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night