




MET2B rabbit pAb

Catalog No	YP-Ab-09159
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	METTL2B
Immunogen	Synthesized peptide derived from human MET2B AA range: 134-184
Specificity	This antibody detects endogenous levels of MET2B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MET2B;METTL2B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	function:Probable methyltransferase.,similarity:Belongs to the methyltransferase superfamily. METL family.,
Background	This gene is a member of a family of methyltransferases that share homology with, but are distinct from, the UbiE family of methyltransferases. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4° over night