



MYNN rabbit pAb


Catalog No	YP-Ab-09160
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MYNN OSZF ZBTB31 SBBIZ1
Immunogen	Synthesized peptide derived from human MYNN AA range: 390-440
Specificity	This antibody detects endogenous levels of MYNN at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MYNN;myoneurin;OSZF;SBBIZ1;ZBTB31;MYNN OSZF ZBTB31 SBBIZ1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Mainly expressed in the neuromuscular system. Located in and around synaptic myonuclei in adult muscle. Expression is dysregulated after nerve injury. Also found in the testis, ovary and placenta.
Function	similarity:Belongs to the krueppel C2H2-type zinc-finger protein family., similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 8 C2H2-type zinc fingers.,tissue specificity:Mainly expressed in the neuromuscular system. Located in and around synaptic myonuclei in adult muscle. Expression is dysregulated after nerve injury. Also found in the testis, ovary and placenta.,
Background	This gene encodes a member of the BTB/POZ and zinc finger domain-containing protein family that are involved in the control of gene expression. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 14. [provided by RefSeq, Jun 2010],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4°over night