




CJ076 rabbit pAb

Catalog No	YP-Ab-09161
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C10orf76
Immunogen	Synthesized peptide derived from human CJ076 AA range: 135-185
Specificity	This antibody detects endogenous levels of CJ076 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CJ076;C10orf76
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass membrane protein . Cytoplasm . The majority of ARMD3 is cytosolic, with a portion colocalizing with GBF1 at juxtanuclear Golgi sites. .
Tissue Specificity	
Function	similarity:Belongs to the UPF0668 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night