




F16A1 rabbit pAb

Catalog No	YP-Ab-09162
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM160A1
Immunogen	Synthesized peptide derived from human F16A1 AA range: 182-232
Specificity	This antibody detects endogenous levels of F16A1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F16A1;FAM160A1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	sequence caution:Contaminating sequence. Potential poly-A sequence., sequence caution:Translation N-terminally extended.,similarity:Belongs to the UPF0518 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night