




CP048 rabbit pAb

Catalog No	YP-Ab-09164
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C16orf48 UNQ6410/PRO21183
Immunogen	Synthesized peptide derived from human CP048 AA range: 7-57
Specificity	This antibody detects endogenous levels of CP048 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CP048;C16orf48 UNQ6410/PRO21183
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasmic microtubule,microtubule cytoskeleton,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night