




T150C rabbit pAb

Catalog No	YP-Ab-09166
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM150C
Immunogen	Synthesized peptide derived from human T150C AA range: 36-86
Specificity	This antibody detects endogenous levels of T150C at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T150C; TMEM150C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein . Localizes at the plasma membrane. A portion co-localizes with LAMP1 lysosomal marker (PubMed:25608530). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night