



TM41A rabbit pAb

Catalog No	YP-Ab-09178
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMEM41A UNQ168/PRO194
Immunogen	Synthesized peptide derived from human TM41A AA range: 11-61
Specificity	This antibody detects endogenous levels of human TM41A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM41A;TMEM41A UNQ168/PRO194;Transmembrane protein 41A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

