



CB066 rabbit pAb

Catalog No	YP-Ab-09182
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C2orf66 UNQ6411/PRO21186
Immunogen	Synthesized peptide derived from human CB066 AA range: 17-67
Specificity	This antibody detects endogenous levels of human CB066
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CB066;C2orf66 UNQ6411/PRO21186;Uncharacterized protein C2orf66
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

