



TPC rabbit pAb

Catalog No	YP-Ab-09187
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SLC25A19 DNC MUP1
Immunogen	Synthesized peptide derived from human TPC AA range: 74-124
Specificity	This antibody detects endogenous levels of human TPC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Papillary thyroid carcinoma-encoded protein;Coiled-coil domain-containing protein 6;coiled coil domain containing 6;TPC;SLC25A19 DNC MUP1; Mitochondrial thiamine pyrophosphate carrier (Mitochondrial uncoupling protein 1) (Solute carrier family 25 member 19)
Observed Band	65 kDa
Calculated Molecular Weight	474 aa, 53 kDa
Cell Pathway	Mitochondrion inner membrane; Multi-pass membrane protein.
Tissue Specificity	Expressed in all tissues examined except for placenta. Highest levels in colon, kidney, lung, testis, spleen, and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

