



YN4349

CNTD1 rabbit pAb

Catalog No	YP-Ab-09189
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CNTD1 CNTD
Immunogen	Synthesized peptide derived from human CNTD1 AA range: 278-328
Specificity	This antibody detects endogenous levels of human CNTD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNTD;CNTD1;CNTD1 CNTD;Cyclin N-terminal domain-containing protein 1
Observed Band	35-37 kDa
Calculated Molecular Weight	330 aa, 37 kDa
Cell Pathway	Nucleus . Cytoplasm . Chromosome . Shuttles between the nucleus and cytoplasm in a stage-specific manner of prophase I cells. Co-localized at crossover sites with PRR19. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images