



MLAS1 rabbit pAb

Catalog No	YP-Ab-09191
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MLLT4-AS1 C6orf124
Immunogen	Synthesized peptide derived from human MLAS1 AA range: 11-61
Specificity	This antibody detects endogenous levels of human MLAS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MLAS1;MLLT4-AS1 C6orf124;Putative uncharacterized protein MLLT4-AS1 (MLLT4 antisense RNA 1) (MLLT4 antisense gene protein 1) (Protein HGC6.4)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

