



EPMIP rabbit pAb

Catalog No	YP-Ab-09192
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	EPM2AIP1 KIAA0766 My007
Immunogen	Synthesized peptide derived from human EPMIP AA range: 180-230
Specificity	This antibody detects endogenous levels of human EPMIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EPM2AIP 1;EPM2A (laforin) interacting protein 1;EPM2A (laforin) IP1;EPM2A-interacting protein 1;Epm2aip1;EPMIP_HUMAN;FLJ1127;KIAA766;Laforin interacting protein 1;Laforin-interacting protein;My7;EPMIP;EPM2AIP1 KIAA0766 My007;EPM2A-interacting protein 1 (Laforin-interacting protein)
Observed Band	
Calculated Molecular Weight	7kD
Cell Pathway	Endoplasmic reticulum .
Tissue Specificity	Expressed in heart, brain, placenta, liver, pancreas, kidney and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images