



WBP4 rabbit pAb

Catalog No	YP-Ab-09193
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WBP4 FBP21 FNBP21
Immunogen	Synthesized peptide derived from human WBP4 AA range: 235-285
Specificity	This antibody detects endogenous levels of human WBP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FBP21;FNBP21;Formin-binding protein 21;WBP-4;WW domain-binding protein 4;WBP4;WBP4 FBP21 FNBP21;WW domain-binding protein 4 (WBP-4) (Formin-binding protein 21) (WW domain-containing-binding protein 4)
Observed Band	
Calculated Molecular Weight	5-55 kDa
Cell Pathway	Nucleus . Nucleus speckle .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images