



SRPRB rabbit pAb

Catalog No	YP-Ab-09194
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SRPRB PSEC0230
Immunogen	Synthesized peptide derived from human SRPRB AA range: 174-224
Specificity	This antibody detects endogenous levels of human SRPRB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	APMCF1;Protein APMCF1;SR beta;SRPRB;SRPRB PSEC0230;Signal recognition particle receptor subunit beta (SR-beta) (Protein APMCF1)
Observed Band	30 kDa
Calculated Molecular Weight	30 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

