



FUT11 rabbit pAb

Catalog No	YP-Ab-09195
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FUT11
Immunogen	Synthesized peptide derived from human FUT11 AA range: 131-181
Specificity	This antibody detects endogenous levels of human FUT11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Alpha-(1;3)-fucosyltransferase 11;EC:2.4.1.-;Fuc TXI;FucT XI;FucTXI;FUT11;3)-fucosyltransferase 11 (EC 2.4.1.-) (Fucosyltransferase XI) (Fuc-TXI) (FucT-XI) (Galactoside 3-L-fucosyltransferase 11) (Fucosyltransferase 11)
Observed Band	56 kDa
Calculated Molecular Weight	492 aa, 56 kDa
Cell Pathway	Golgi apparatus, Golgi stack membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

