



1A66 rabbit pAb

Catalog No	YP-Ab-09198
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HLA-A HLAA
Immunogen	Synthesized peptide derived from human 1A66 AA range: 22-72
Specificity	This antibody detects endogenous levels of human 1A66
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1A66;HLA-A HLAA;HLA class I histocompatibility antigen;A-66 alpha chain (Aw-66) (HLA class I histocompatibility antigen;A-10 alpha chain) (MHC class I antigen A*66)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

