



TNPO3 rabbit pAb

Catalog No	YP-Ab-09200
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TNPO3 IPO12
Immunogen	Synthesized peptide derived from human TNPO3 AA range: 601-651
Specificity	This antibody detects endogenous levels of human TNPO3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Imp12;Importin 12;Importin-12;IPO12;transportin 3;TNPO3;TNPO3 IPO12; Transportin-3 (Importin-12) (Imp12) (Transportin-SR) (TRN-SR)
Observed Band	104 kDa
Calculated Molecular Weight	975 aa, 110 kDa
Cell Pathway	Nucleus envelope . Cytoplasm . Localizes to the nuclear envelope and annulate lamellae, which consists in stacks of endoplasmic reticulum membranes containing a high density of nuclear pores. .
Tissue Specificity	Expressed in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

