



YN4361

UN93A rabbit pAb

Catalog No	YP-Ab-09201
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UNC93A
Immunogen	Synthesized peptide derived from human UN93A AA range: 312-362
Specificity	This antibody detects endogenous levels of human UN93A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UN93A;UNC93A;Protein unc-93 homolog A (HmUnc-93A) (Unc-93A)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in testis, small intestine, spleen, prostate and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images