



NUD15 rabbit pAb

Catalog No	YP-Ab-09202
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NUDT15 MTH2
Immunogen	Synthesized peptide derived from human NUD15 AA range: 62-112
Specificity	This antibody detects endogenous levels of human NUD15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NUD15;NUDT15 MTH2;Probable 8-oxo-dGTP diphosphatase NUDT15 (8-oxo-dGTPase NUDT15) (EC 3.6.1.55) (7;8-dihydro-8-oxoguanine-triphosphatase NUDT15) (MutT homolog 2) (MTH2) (Nucleoside diphosphate-linked moiety X motif 15)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

