



PRRT1 rabbit pAb

Catalog No	YP-Ab-09204
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRRT1 C6orf31 NG5
Immunogen	Synthesized peptide derived from human PRRT1 AA range: 19-69
Specificity	This antibody detects endogenous levels of human PRRT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRRT1;Proline-rich transmembrane protein 1;NG5;DSPD1;Dispanin subfamily D member 1;PRRT1 C6orf31 NG5;Proline-rich transmembrane protein 1 (Dispanin subfamily D member 1) (DSPD1)
Observed Band	35-40 kDa
Calculated Molecular Weight	225aa,23 kDa; 306aa,31 kDa
Cell Pathway	Cell membrane ; Single-pass type II membrane protein . Cell junction, synapse .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

