



XAF1 rabbit pAb

Catalog No	YP-Ab-09207
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	XAF1 BIRC4BP XIAPAF1
Immunogen	Synthesized peptide derived from human XAF1 AA range: 47-97
Specificity	This antibody detects endogenous levels of human XAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BIRC4 binding protein;BIRC4-binding protein;BIRC4BP;XIAP associated factor 1;XIAPAF1;XAF1;XAF1 BIRC4BP XIAPAF1;XIAP-associated factor 1 (BIRC4-binding protein)
Observed Band	35 kDa
Calculated Molecular Weight	301 aa, 35 kDa
Cell Pathway	Cytoplasm. Nucleus. Mitochondrion. Found in the cytoplasm and nucleus of placental syncytiotrophoblasts. Translocates to mitochondria upon TNF-alpha treatment.; [Isoform 1]: Nucleus.; [Isoform 5]: Nucleus.
Tissue Specificity	Widely expressed. Expression is frequently down-regulated in cancer cell lines. Isoform 5 is widely expressed. Expressed in placenta (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

