



IGDC4 rabbit pAb

Catalog No	YP-Ab-09208
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	IGDCC4 DDM36 KIAA1628 NOPE
Immunogen	Synthesized peptide derived from human IGDC4 AA range: 56-106
Specificity	This antibody detects endogenous levels of human IGDC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IGDC4;IGDCC4 DDM36 KIAA1628 NOPE;Immunoglobulin superfamily DCC subclass member 4 (Neighbor of punc e11) (Protein DDM36) (hDDM36)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

