



FBSP1 rabbit pAb

Catalog No	YP-Ab-09210
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXO45 FBX45
Immunogen	Synthesized peptide derived from human FBSP1 AA range: 120-170
Specificity	This antibody detects endogenous levels of human FBSP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FBSP1;FBXO45 FBX45;F-box/SPRY domain-containing protein 1 (F-box only protein 45) (hFbxo45)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane . Cell junction, synapse, presynaptic cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

