



SC5A7 rabbit pAb

Catalog No	YP-Ab-09213
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC5A7 CHT1
Immunogen	Synthesized peptide derived from human SC5A7 AA range: 149-199
Specificity	This antibody detects endogenous levels of human SC5A7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SC5A7;SLC5A7 CHT1;High affinity choline transporter 1 (Hemicholinium-3-sensitive choline transporter) (CHT) (Solute carrier family 5 member 7)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . Cell membrane . Cell junction, synapse . Localized at the neuromuscular junction. .
Tissue Specificity	Expressed in putamen, spinal cord and medulla. Specific for cholinergic neurons.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

