



NHP2 rabbit pAb

Catalog No	YP-Ab-09217
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NHP2 NOLA2 HSPC286
Immunogen	Synthesized peptide derived from human NHP2 AA range: 70-120
Specificity	This antibody detects endogenous levels of human NHP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NOLA2;NHP2P;NHP 2;HSPC286;NHP2;NHP2 NOLA2 HSPC286;H/ACA ribonucleoprotein complex subunit 2 (Nucleolar protein family A member 2) (snoRNP protein NHP2)
Observed Band	17-20 kDa
Calculated Molecular Weight	17 kDa
Cell Pathway	Nucleus, nucleolus. Nucleus, Cajal body. Also localized to Cajal bodies (coiled bodies).
Tissue Specificity	Expressed in brain, colon, heart, kidney, ovary, pancreas, placenta, prostate, skeletal muscle, small intestine, spleen, testis and thymus. Also expressed at lower levels in the liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images