



ZN410 rabbit pAb

Catalog No	YP-Ab-09218
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZNF410 APA1
Immunogen	Synthesized peptide derived from human ZN410 AA range: 113-163
Specificity	This antibody detects endogenous levels of human ZN410
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Another partner for ARF 1;APA 1;APA1;Clones 23667 and 23775 zinc finger protein;D12Ertd748e;Zfp41;Zinc finger protein 41;Zinc finger protein APA 1;Zinc finger protein APA-1;ZN41_HUMAN;Znf41;ZN410;ZNF410 APA1;Zinc finger protein 410 (Another partner for ARF 1) (Zinc finger protein APA-1)
Observed Band	
Calculated Molecular Weight	52kD
Cell Pathway	Nucleus . Chromosome . Directly binds to the sequence motif 5'-CATCCCATAATA-3' .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

