



TRUB1 rabbit pAb

Catalog No	YP-Ab-09220
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRUB1 PUS4
Immunogen	Synthesized peptide derived from human TRUB1 AA range: 116-166
Specificity	This antibody detects endogenous levels of human TRUB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:5.4.99.-;EC:5.4.99.25;Pseudouridylate synthase TRUB1;Psi55 synthase TRUB1;PUS4;TRUB1;TRUB1 PUS4;Probable tRNA pseudouridine synthase 1 (EC 5.4.99.-)
Observed Band	37 kDa
Calculated Molecular Weight	349 aa, 37 kDa
Cell Pathway	Nucleus . Cytoplasm, cytosol . Catalyzes pseudouridylation of mRNAs in the nucleus. .
Tissue Specificity	Highly expressed in heart, skeletal muscle and liver. Expressed at lower levels in lung, small intestine, kidney and spleen.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

