



TMIE rabbit pAb

Catalog No	YP-Ab-09221
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMIE
Immunogen	Synthesized peptide derived from human TMIE AA range: 12-62
Specificity	This antibody detects endogenous levels of human TMIE
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	deafness autosomal recessive 6;DFNB6;transmembrane inner ear; transmembrane inner ear expressed protein;transmembrane inner ear protein; TMIE
Observed Band	
Calculated Molecular Weight	17kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in many tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

