



R113A rabbit pAb

Catalog No	YP-Ab-09223
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RNF113A RNF113 ZNF183
Immunogen	Synthesized peptide derived from human R113A AA range: 76-126
Specificity	This antibody detects endogenous levels of human R113A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	R113A;RNF113A RNF113 ZNF183;RING finger protein 113A (Zinc finger protein 183)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nucleus speckle . Colocalizes with ASCC2 in nuclear foci after DNA damage by alkylating agents. In the absence of DNA damage, colocalizes with the spliceosome components SNRNP200/BRR2 and PRPF8 in nuclear speckles. .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

