



ZBED4 rabbit pAb

Catalog No	YP-Ab-09224
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZBED4 KIAA0637
Immunogen	Synthesized peptide derived from human ZBED4 AA range: 10-60
Specificity	This antibody detects endogenous levels of human ZBED4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZBED4;ZBED4 KIAA0637;Zinc finger BED domain-containing protein 4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm . Photoreceptor inner segment .
Tissue Specificity	Expressed in testis, heart, lung, and weakly expressed in brain, liver, muscle, placenta and small intestine (PubMed:19369242). Expressed in the retina, found in the cone photoreceptors, Mueller cells, cone pedicles and in the innermost retinal layer (PubMed:19369242).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images