



FXL20 rabbit pAb

Catalog No	YP-Ab-09225
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXL20 FBL2
Immunogen	Synthesized peptide derived from human FXL20 AA range: 348-398
Specificity	This antibody detects endogenous levels of human FXL20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FXL20;FBXL20 FBL2;F-box/LRR-repeat protein 20 (F-box and leucine-rich repeat protein 20) (F-box/LRR-repeat protein 2-like)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

