



TMF1 rabbit pAb

Catalog No	YP-Ab-09227
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMF1 ARA160
Immunogen	Synthesized peptide derived from human TMF1 AA range: 498-548
Specificity	This antibody detects endogenous levels of human TMF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMF1;TMF;TATA element modulatory factor;ARA160;Androgen receptor-associated protein of 160 kDa;TMF1 ARA160;TATA element modulatory factor (TMF) (Androgen receptor coactivator 160 kDa protein) (Androgen receptor-associated protein of 160 kDa)
Observed Band	150-160 kDa
Calculated Molecular Weight	123 kDa
Cell Pathway	Cytoplasm. Nucleus. Golgi apparatus membrane. Concentrated at the budding structures localized at the tips of cisternae.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

