



BRPF1 rabbit pAb

Catalog No	YP-Ab-09228
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	BRPF1 BR140
Immunogen	Synthesized peptide derived from human BRPF1 AA range: 527-577
Specificity	This antibody detects endogenous levels of human BRPF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BR140;Bromodomain and PHD finger-containing protein 1;Peregrin;Protein Br140;BRPF1;BRPF1 BR140;Peregrin (Bromodomain and PHD finger-containing protein 1) (Protein Br140)
Observed Band	140 kDa
Calculated Molecular Weight	1220 aa, 138 kDa
Cell Pathway	Nucleus . Chromosome . Cytoplasm . Localization to the nucleus depends on KAT6A, ING5 and MEAF6 (PubMed:18794358, PubMed:27939640). Localizes to transcription start sites (PubMed:24065767). .
Tissue Specificity	High levels in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

