



YN4389

UB2J1 rabbit pAb

Catalog No	YP-Ab-09229
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UBE2J1 NCUBE1 CGI-76 HSPC153 HSPC205
Immunogen	Synthesized peptide derived from human UB2J1 AA range: 39-89
Specificity	This antibody detects endogenous levels of human UB2J1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UB2J1;UBE2J1 NCUBE1 CGI-76 HSPC153 HSPC205;Ubiquitin-conjugating enzyme E2 J1 (EC 6.3.2.19) (Non-canonical ubiquitin-conjugating enzyme 1) (NCUBE-1) (Yeast ubiquitin-conjugating enzyme UBC6 homolog E) (HsUBC6e)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images