



ZDHC4 rabbit pAb

Catalog No	YP-Ab-09230
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZDHHC4 ZNF374 DC1 UNQ5787/PRO19576
Immunogen	Synthesized peptide derived from human ZDHC4 AA range: 109-159
Specificity	This antibody detects endogenous levels of human ZDHC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZDHC4;ZDHHC4 ZNF374 DC1 UNQ5787/PRO19576;Probable palmitoyltransferase ZDHHC4 (EC 2.3.1.-) (Zinc finger DHHC domain-containing protein 4) (DHHC-4) (Zinc finger protein 374)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

